




Probe	CD274/CEN9q Dual Color FISH Probe(V.2)
Product No.	GC197
CD274	
Probe Size	Approximately 150kb
Fluorophore	TexRed
Location	9p24
Also known as	B7-H; B7H1; PDL1; PD-L1; PDCD1L1; PDCD1LG1

CEN9q	
Probe Size	Approximately 470kb
Fluorophore	FITC
Location	9q21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.