

Probe DPYD/CEN1p Dual Color FISH Probe
Product No. GC192


DPYD

Probe Size Approximately 800kb
Fluorophore TexRed
Location 1p22
Also known as DHP; DPD; DHPDHASE




CEN1p

Probe Size Approximately 780kb
Fluorophore FITC
Location 1p13.3

DPYD



CEN1p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.