

Probe EGR1/CEN5q Dual Color FISH Probe
Product No. GC186


EGR1

Probe Size Approximately 260kb
Fluorophore TexRed
Location 5q31.1
Also known as TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268

CEN5q

Probe Size Approximately 700kb
Fluorophore FITC
Location 5q11.2

EGR1





~260kb

CEN5q



~700kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.