



Probe CD274/CEN9q Dual Color FISH Probe
Product No. GC141

CD274

Probe Size Approximately 230kb
Fluorophore TexRed
Location 9p24
Also known as B7-H; B7H1; PDL1; PD-L1; PDCD1L1; PDCD1LG1

CEN9q

Probe Size Approximately 470kb
Fluorophore FITC
Location 9q21



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.

Reference

Mod Pathol. 2016 Dec 9. doi: 10.1038/modpathol.2016.202. [Epub ahead of print]

Overexpression and gene amplification of PD-L1 in cancer cells and PD-L1+ immune cells in Epstein-Barr virus-associated gastric cancer: the prognostic implications