

Probe FANCA/CEN16q Dual Color FISH Probe
Product No. GC205

FANCA

Probe Size Approximately 360kb
Fluorophore TexRed
Location 16q24.3
Also known as FA; FA1; FAA; FAH; FA-H; FACA; FANCH

CEN16q

Probe Size Approximately 700kb
Fluorophore FITC
Location 16q12.1

FANCA

89.2 89.4 89.6 89.8 90.0 90.2 90.4 (MB)





~360kb

CEN16q

48.6 48.8 49.0 49.2 49.4 49.6 49.8 (MB)



~700kb



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.