

Probe CDKN1B/CEN12p Dual Color FISH Probe
Product No. GC209


CDKN1B

Probe Size Approximately 140kb
Fluorophore TexRed
Location 12p13.1
Also known as KIP1; MEN4; CDKN4; MEN1B; P27KIP1




CEN12p

Probe Size Approximately 530kb
Fluorophore FITC
Location 12p11.22

CDKN1B



CEN12p



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.