




Probe c-met/CEN7 Dual Color FISH Probe
Product No. GC215

c-met
Probe Size Approximately 270kb
Fluorophore TexRed
Location 7q31
Also known as HGFR; AUTS9; RCCP2; MET; DFN97

CEN7(D7Z1)
Probe Size Approximately 3Mb
Fluorophore FITC
Location 7p11.1-q11.1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.