

Probe TYMP/CEN22q Dual Color FISH Probe
Product No. GC153

TYMP

Probe Size Approximately 150kb

Fluorophore TexRed

Location 22q13

Also known as TP; ECGF; ECGF1; MNGIE; MEDPS1; MTDPS1; PDECGF; hPD-ECGF


CEN22q

Probe Size Approximately 620kb




Fluorophore FITC

Location 22q11.22

TYMP



CEN22q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.