






Probe HER2/TOP1/CEN17 TriColor FISH Probe
Product No. GC248

HER2
Probe Size Approximately 510kb
Fluorophore FITC
Location 17q12
Also known as NEU; NGL; ERBB2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2)

TOP1
Probe Size Approximately 420kb
Fluorophore TAMRA (Orange)
Location 20q12

CEN17 (D17Z1)
Fluorophore DEAC
Location 17p11.1-17q11.1





representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.