



Probe	BRAF1 Split Dual Color FISH Probe
Product No.	SP060
BRAF1	
Fluorophore	TexRed, FITC
Location	7q34
Also known as	NS7; B-raf; BRAF1; RAFB1; B-RAF1
Storage	4°C



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.