

Probe PD-L1/PD-L2 Split FISH Probe
Product No SP075




PD-L1/PDL2

Fluorophore TexRed、FITC

Location 9p24

Also known as B7-H; B7H1; PDL1; CD274; PDCD1L1; PDCD1LG1/ B7DC; Btdc; PDL2; CD273; PDCD1LG2; PDCD1L2; bA574F11.2

Storage 4°C



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.