




Probe VHL Split Dual Color FISH Probe
Product No. SP090

VHL
Fluorophore TexRed, FITC
Location 3p25.3
Also known as RCA1; VHL1; pVHL; HRCA1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.