



Probe	IRF4, DUSP22 Split Dual Color FISH Probe
Product No.	SP079
IRF4, DUSP22	
Fluorophore	TexRed、FITC
Location	6p25.3
Also known as	MUM1; LSIRF; SHEP8; NF-EM5 VHX; JKAP; JSP1; MKPX; JSP-1; MKP-x; LMWDSP2; LMW-DSP2



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.