

Probe IGH/FOXP1, DY Translocation Dual Color FISH Probe
Product No. TR017


IGH

Probe Size Approximately 1800kb
Fluorophore FITC
Location 14q32
Also known as IGD1; IGH@; IGHJ; IGHV; IGHD@; IGHJ@; IGHV@; IGH.1@; IGHDY1




FOXP1

Probe Size Approximately 1180kb
Fluorophore TexRed
Location 3p14.1
Also known as MFH; QRF1; 12CC4; hFKH1B; HSPC215

IGH



FOXP1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.