

Probe MYB/NFIB, DY Translocation Dual Color FISH Probe
Product No. TR031


MYB

Probe Size Approximately 1010kb
Fluorophore TexRed
Location 6q22-23
Also known as efg; Cmyb; c-myb; c-myb_CDS




NFIB

Probe Size Approximately 1090kb
Fluorophore TexRed
Location 9p24.1
Also known as CTF; NF1-B; NFI-B; NFIB2; NFIB3; NF-I/B; NFI-RED; HMGIC/NFIB

MYB



NFIB



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.