

Probe CSF1/COL6A3, DY Translocation FISH Probe
Product No. TR032


CSF1

Probe Size Approximately 850kb
Fluorophore FITC
Location 1p13.3
Also known as MCSF; CSF-1




COL6A3

Probe Size Approximately 680kb
Fluorophore TexRed
Location 2q37.3
Also known as DYT27; UCMD1; BTHLM1

CSF1



COL6A3



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.