

Probe CSF1/NBPF9, SY Translocation FISH Probe
Product No. TR042


CSF1

Probe Size Approximately 660kb
Fluorophore FITC
Location 1p13.3
Also known as MCSF; CSF-1




NBPF9

Probe Size Approximately 350kb
Fluorophore TexRed
Location 1q21.2
Also known as AE01

CSF1



NBPF9



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.