

Probe RREB1/MKL2, SY Translocation FISH Probe
Product No. TR043


RREB1

Probe Size Approximately 610kb
Fluorophore FITC
Location 6p24.3
Also known as HNT; FINB; LZ321; Zep-1; RREB-1




MKL2

Probe Size Approximately 560kb
Fluorophore TexRed
Location 16p13.12
Also known as MRTFB; MRTF-B; NPD001

RREB1



MKL2



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.