

Probe NUT/BRD4, SY Translocation FISH Probe
Product No. TR044


NUT

Probe Size Approximately 530 kb
Fluorophore TexRed
Location 15q14
Also known as NUTM1; FAM22H; C15orf55




BRD4

Probe Size Approximately 730kb
Fluorophore FITC
Location 19p13.1
Also known as CAP; MCAP; HUNK1; HUNKI

NUT



BRD4



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.