

Probe TFEB/MALAT1, SY Dual Color Translocation FISH Probe
Product No. TR046


TFEB

Probe Size Approximately 470 kb
Fluorophore FITC
Location 6p21.2
Also known as TCFEB; BHLHE35; ALPHATFEB




MALAT1

Probe Size Approximately 640 kb
Fluorophore TexRed
Location 11q13.1
Also known as HCN; NEAT2; PRO2853; LINC00047; NCRNA00047

TFEB



MALAT1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.