

Probe KIF5B/RET, SY Translocation Dual Color FISH Probe(V.2)
Product No. TR050


KIF5B

Probe Size Approximately 720kb
Fluorophore TexRed
Location 10p11.2
Also known as KNS; KINH; KNS1; UKHC; HEL-S-61




RET

Probe Size Approximately 630kb
Fluorophore FITC
Location 10q11.2
Also known as PTC; MTC1; HSCR1; MEN2A; MEN2B; CDHF12; CDHR16; RET-ELE1

KIF5B



RET



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.