





Probe MECOM/RPN1, DY Translocation Dual Color FISH Probe
Product No. TR055

MECOM

Probe Size Approximately 1090kb
Fluorophore TexRed
Location 3q26.2
Also known as EVI1; MDS1; KMT8E; PRDM3; RUSAT2; MDS1-EVI1; AML1-EVI-1

RPN1

Probe Size Approximately 1350kb
Fluorophore FITC
Location 3q21.3
Also known as OST1; RBPH1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.