

Probe TCF3/PBX1, SY Translocation Dual Color FISH Probe
Product No. TR056


TCF3

Probe Size Approximately 840kb
Fluorophore FITC
Location 19p13.3
Also known as E2A; E47; p75; AGM8; ITF1; VDIR; AGM8A; AGM8B; TCF-3; bHLHb21




PBX1

Probe Size Approximately 610kb
Fluorophore TexRed
Location 1q23.3
Also known as CAKUHD

TCF3



PBX1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.