




Probe API2/MALT1, DY Translocation Dual Color FISH Probe
Product No. TR057

API2

Probe Size Approximately 1350kb
Fluorophore FITC
Location 11q21
Also known as AIP1; BIRC3; MIHC; CIAP2; HAIP1; HIAP1; IAP-1; MALT2; RNF49; c-IAP2

MALT1

Probe Size Approximately 890kb
Fluorophore TexRed
Location 18q21
Also known as MLT; MLT1; IMD12; PCASP1



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.