

Probe PML/RARA,DY Translocation Dual Color FISH Probe
Product No. TR007


RARA

Probe Size Approximately 1200kb
Fluorophore FITC
Location 17q21.2
Also known as RAR; NR1B1




PML

Probe Size Approximately 530kb
Fluorophore TexRed
Location 15q24.1
Also known as MYL; RNF71; PP8675; TRIM19

RARA



PML



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.