

Probe VCL/ALK, DY Translocation Dual Color FISH Probe
Product No. TR021


VCL

Probe Size Approximately 1140kb
Fluorophore TexRed
Location 10q22.2
Also known as MV; MVCL; CMD1W; CMH15; HEL114




ALK

Probe Size Approximately 1220kb
Fluorophore FITC
Location 2p23
Also known as CD246; NBLST3

VCL



ALK



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.