

Probe MST1R/CEN3q Dual Color FISH Probe  
Product No. GC098



**MST1R**

Probe Size Approximately 150kb  
Fluorophore TexRed  
Location 3p21.3  
Also known as RON; PTK8; CD136; CDw136




**CEN3q**

Probe Size Approximately 500kb  
Fluorophore FITC  
Location 3q12.1

MST1R



CEN3q



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.