




Probe HMGA2/NCOR2, SY Translocation Dual Color FISH Probe  
Product No. TR052

**HMGA2**

Probe Size Approximately 650kb  
Fluorophore TexRed  
Location 12q14.3  
Also known as BABL; LIPO; SRS5; HMGIC; HMGI-C; STQTL9

**NCOR2**

Probe Size Approximately 660kb  
Fluorophore FITC  
Location 12q24.31  
Also known as SRC2; TIF2; GRIP1; KAT13C; NCoA-2; bHLHe75



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.