





Probe THBS1/ADGRF5, SY Translocation Dual Color FISH Probe
Product No. TR051

THBS1
Probe Size Approximately 700kb
Fluorophore FITC
Location 15q14
Also known as TSP; THBS; TSP1; TSP-1; THBS-1

ADGRF5
Probe Size Approximately 670kb
Fluorophore TexRed
Location 6p12.3
Also known as GPR116; KPG_001



Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.