
Probe	EML4/ALK, DY Translocation Dual Color FISH Probe
Product No.	TR001


EML4

Probe Size	Approximately 1130kb
Fluorophore	FITC
Location	2p22-p21
Also known as	C2orf2; ELP120; EMAP-4; EMAPL4; ROPP120

ALK

Probe Size	Approximately 1220kb
Fluorophore	TexRed
Location	2p23
Also known as	CD246; NBLST3





Representative images of normal human cell (lymphocyte) stain with the dual color FISH probe. The left image is chromosomes at metaphase, and the right image is an interphase nucleus.